Summary and Explanation of the Device

Standard methods of biological specimen storage and transport have several disadvantages, such as low-temperature or low-volume requirements. ViveST™ provides the ability to store up to 1mL of viable biological specimen at ambient temperature for use in specimen storage and transportation.^{1,2} This device utilizes an absorbent matrix on which a biological specimen can be loaded and dried. The matrix is housed in a screw-top tube such that the sample is self-contained during storage and shipping. The tube also contains a DriCap® Desiccator that acts as a quality control indicator for dryness during shipping and storage.

Device Principle

Using the ViveST Sample Storage and Transportation Device, plasma is added and allowed to dry on an absorbent matrix. During the drying process, water in the plasma is evaporated but biological analytes including proteins, virions and nucleic acids remain in the matrix. The samples can then be stored or transported as needed in the dried state. The samples are then reconstituted into a liquid state such that the dried analytes are released from the matrix and contained in molecular-grade water. Reconstituted samples can then be processed to determine HIV viral load or HIV genotype.

Intended Use/Indications for Use

This product is intended for drying and reconstitution of human plasma samples containing the HIV virus collected with K₂EDTA. The reconstituted samples are intended for processing using the Siemens VERSANT® HIV-1 RNA 3.0 Assay (bDNA) or TRUGENE® HIV-1 Genotyping Assay.

Materials Provided

Store all components at Room Temperature (15-25°C, out of direct sunlight).

Component	Description	Quantity
ViveST	Tube containing fibrous matrix and DriCap Desiccator	24


Materials Required, but Not Provided

Item	Description
Pipette	1000 µL
Pipette Tips	1000 µL, sterile DNase- and RNase-free, with aerosol barrier
Water	Molecular-grade
Sterile Syringes	BD 3 mL Luer-Lock™, individually-wrapped
Tubes	At least 1 mL with a screw-cap, for sample recovery/storage
Rack	For tubes used in recovery

Workstation	Biological safety cabinet or laminar flow hood
Personal Protection	Powder-free disposable gloves, laboratory coat, safety glasses
Waste	Biohazard waste container
Hygrometer	Measures laboratory environmental humidity

Warnings and Precautions

- For *in vitro* Diagnostic Use.
- For Export Only.
- This kit supplies sufficient consumables and reagents for the number of tests indicated.
- Do not use the components of this kit beyond the expiration date printed on the box.
- Use only a ViveST containing a DriCap Desiccator with a deep blue color. Do not use a ViveST if the DriCap Desiccator appears light-blue, white or pink in color.
- Use good laboratory practices and universal precautions relating to the prevention of transmission of blood borne pathogens.^{3,4}
 - Use disposable powder-free gloves to handle all materials as though capable of transmitting infectious agents.
 - Immediately clean any spills suspected of potentially containing infectious agents with 0.5% w/v sodium hypochlorite (10% v/v bleach).
 - Dispose of all specimens and materials that come in contact with specimens as though they contain infectious agents.
 - In the event that materials known or suspected of containing infectious agents are ingested or come in contact with open lacerations, lesions or mucous membranes (eyes, nasal passages, etc.), consult with a physician immediately.

Symbology Table			
Symbol	Meaning	Symbol	Meaning
15°C  25°C	Temperature limitation		Contains sufficient for n tests
	Use by		Catalogue number
	Batch code		Consult instructions for use
	Do not reuse		Caution/Attention, see instructions for use
	Manufacturer		Date of Manufacture
	European Conformity		For <i>In Vitro</i> Diagnostic Use

Sample Loading Procedure

1. Place desired number of ViveST into a biological safety cabinet or laminar flow hood.
2. Observe the DriCap Desiccator in each ViveST. Use only a ViveST containing a DriCap Desiccator that is deep blue in color.

Warning - Do not use if the DriCap Desiccator appears light-blue, white or pink in color.

3. Label the **cap** of each ViveST with a sample identifier.

- The matrix is attached to the cap; therefore, the sample identifier should be placed on the cap. Labeling only the tube may lead to sample confusion.
4. Remove the cap from each ViveST and place on the surface of the hood such that the absorbent matrix is pointing upward.
 - Retain the ViveST tube in the hood out of the way.
 5. Slowly pipette 1 mL of the sample to the top of the ViveST matrix.
 - Initially adding 1-2 drops of sample may help to begin proper absorption.
 - It is acceptable to lightly touch the matrix with the pipette tip.
 6. Air dry the loaded matrix.
 - Do NOT place the matrix into the ViveST tube until the matrix has completely dried.
 - Samples will dry faster if placed in an air ventilation instrument, such as a laminar flow hood or biological safety cabinet, near the circulating vent.
 - Generally, 8 hours are required for drying when using a laminar flow hood. This time may vary due to various environmental factors such as airflow speeds, temperature and humidity. At relative humidity >60%, drying time may need to be extended.
 - For individual laboratory validation, follow the ViveST Drying Time Validation Protocol. For consistent drying results, use a hygrometer to monitor the humidity of the drying area.
 7. Screw the cap onto the ViveST such that the matrix is contained within the tube.

Warning – The matrix must be completely dry before replacing it into the ViveST.
 8. The ViveST containing dried sample may be stored for up to 2 months or shipped at ambient temperature.
 - Ship the ViveST containing the dried sample following 2009 UN Model Regulations and WHO 2009-2010 Guidance.^{5,6}
 9. Refer to the below Sample Recovery Procedure for instructions on reconstitution of the sample.

Sample Recovery Procedure

1. Place the ViveSTs containing dried samples into a clean biological hood.
2. Observe the DriCap Desiccator inside each ViveST. Continue recovery only for a ViveST containing a DriCap Desiccator that is a deep blue color.

Warning - Do not use if the DriCap Desiccator appears light-blue, white or pink in color.
3. Inside the hood, uncap a tube for recovery of each sample and place into the tube rack.
 - Place any recovery tube screw caps in a clean area of the hood until ready for further use. Do NOT discard the caps.
4. Prepare a sterile syringe by removing its wrapping and removing the plunger.
 - Place the removed plungers in a clean area of the hood until ready for further use. Do NOT discard the plungers.
5. Uncap one ViveST and transfer the matrix into the syringe barrel.
 - Pressing the matrix against the inside of the syringe barrel mouth should apply enough pressure to break the matrix free from the cap.
 - Holding the matrix against the syringe barrel with a pipette tip may also aid in releasing the matrix.

6. Place the syringe barrel into a recovery tube.
7. Draw 1 mL of molecular grade water into a pipette tip.
8. Using the pipette tip, gently push the matrix until it is touching the bottom of the syringe barrel. Pipette the water onto the matrix.
 - Do not squeeze the matrix when placing it to the bottom of the syringe barrel or when pipetting the water
 - The water should be completely absorbed by the matrix. If buffer leaks into the recovery tube, aspirate it using the pipette tip and re-pipette onto the matrix.
9. Rest the syringe in the recovery tube and allow the matrix to incubate at room temperature for 10 minutes.
10. Perform steps 4 - 9 for each ViveST.
11. While keeping the tip of each syringe in its recovery tube, insert the syringe plunger into the syringe barrel and depress with firm, even pressure until the plunger has completely compressed the matrix.
 - A volume of approximately 1 mL should be collected in the recovery tube.
12. Remove the syringe from the recovery tube and dispose of all components, including the matrix, in a proper biohazard waste container.
13. Cap and label each recovery tube containing reconstituted sample.
14. Continue with testing or store samples using your standard plasma storage conditions.

Function of the DriCap Desiccator

The DriCap Desiccator is intended for use as quality control of sample dryness as well as the storage and shipping environments. It is NOT intended to aid in drying of the samples on the matrix. The matrix should be completely dried prior to being secured in the ViveST. If the matrix is not properly dried, the DriCap Desiccator will change from a deep blue color to varying shades of light-blue, pink and white. If any of the crystals in the desiccant are light-blue, pink or white, the sample should not be used for further testing due to possible analyte degradation.

Test Limitations

- Validation tests of ViveST have only been performed using human plasma collected with K₂EDTA and processed using the Siemens VERSANT[®] HIV-1 RNA 3.0 Assay (bDNA) and TRUGENE[®] HIV-1 Genotyping Assay.
- The DriCap Desiccator is NOT intended to dry samples. Samples must be completely dry before securing the matrix in the ViveST.

Performance Characteristics

Tests compared plasma samples stored on the ViveST Sample Storage and Transportation Device to the gold standard of frozen plasma for storage and transportation. Studies were conducted using two assays: Siemens VERSANT[®] HIV-1 RNA 3.0 Assay (bDNA) and TRUGENE[®] HIV-1 Genotyping Assay. The specimens selected for testing contained various viral loads, genotypes and drug resistance mutations.

Test	# of Samples
VERSANT	296
TRUGENE	230

Test specimens were selected to meet the standards set forth in the manufacturer's package inserts for each assay. Additionally, the manufacturer's performance data was used as a guideline for the calculation of the

performance data below. For the VERSANT assay, specimens were compared to determine the ability to detect their viral load. For the TRUGENE assay, specimens were compared at the nucleotide level such that the number of concordant nucleic acids (NA) between the ViveST-processed and frozen plasma were determined.

Sensitivity

Sensitivity was calculated as the positive result concordance. For the VERSANT assay, sensitivity was calculated as the number of ViveST-processed specimens that give a detectable result (>50 copies/mL) divided by the number of frozen plasma specimens that give a detectable result.

Test	# ViveST Detectable	# Frozen Detectable	Sensitivity
VERSANT	187	198	94%

For the TRUGENE assay, sensitivity was calculated as the mutation recognition agreement rate, the number of resistance-associated nucleotides concordant between specimens processed on the ViveST and frozen plasma specimens divided by the total number of resistance-associated nucleotides.

Test	Concordant NA	Total NA	Sensitivity
TRUGENE	34,956.0	35,970.0	97%

Specificity

Specificity was calculated as the negative result concordance. For the VERSANT assay, specificity was calculated as the agreement between the number of ViveST-processed specimens that gave an undetectable result (<50 copies/mL) divided by the number of undetectable from frozen plasma.

Test	# ViveST Undetectable	# Frozen Undetectable	Specificity
VERSANT	87	98	89%

For the TRUGENE assay, specificity was calculated as the wild type recognition agreement rate, the number of non-resistance-associated nucleotides concordant between specimens processed on the ViveST and frozen plasma specimens divided by the total number of non-resistance-associated nucleotides.

Test	Concordant NA	Total NA	Specificity
TRUGENE	164,960.5	165,680.0	99%

Accuracy

Accuracy was calculated as the concordance between the paired results for the ViveST-processed plasma and the frozen plasma. Specifically for the VERSANT assay, accuracy was calculated as the percentage of ViveST-processed specimens that showed the same result (detectable or undetectable) as the frozen plasma.

Test	# Samples Concordant	# Samples Tested	Accuracy
VERSANT	274	296	92%

For the TRUGENE assay, accuracy was calculated as the total number of nucleotides concordant between ViveST-processed plasma and frozen plasma divided by the total number of nucleotides assessed.

Test	Concordant NA	Total NA	Accuracy
TRUGENE	199,916.5	201,650.0	99%

References

1. RM Lloyd, Jr., DA Burns, JT Huong, RL Mathis, MA Winters, M Tanner, A De La Rosa, B Yen-Lieberman, W Armstrong, A Taeye, DR McClemon, JL Wetshtein, BM Friedrich, MR Ferguson, W O'Brien, PM Feorino and M Holodniy. Dried-plasma transport using a novel matrix and collection system for human immunodeficiency virus and hepatitis C virus virologic testing. *J Clin Microbiol*, 47(5), May 2009, p. 1491-96.
2. M Zanoni, R Cortes, RS Diaz, MC Sucupira, D Ferreira, LA Inocencio, C Vilhena, C Loveday, RM Lloyd, Jr., M. Holodniy. Comparative effectiveness of dried plasma HIV-1 viral load testing in Brazil using ViveST for sample collection. *J Clin Virol*, 48(4), Dec 2010, p. 245-8.
3. National Committee for Clinical Laboratory Standards. Protection of laboratory workers from infectious disease transmitted by blood, body fluids, and tissue; approved guideline. NCCLS Document M29-A Villanova (PA); NCCLS; 1997 Dec. 90p.
4. Centers for Disease Control and Prevention: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other blood borne pathogens in healthcare settings. *MMWR*, 1988; 37: 377-82, 387-8.
5. United Nations: Recommendations on Transport of Dangerous Goods Model Regulations. UN document, 2009, section 2.6.3.
6. World Health Organization: Guidance on Regulations for the Transport of Infectious Substances 2009-2010. WHO document, 2008.

Contact Information

ViveBio, LLC
1000 Hurricane Shoals Rd
Suite C330
Lawrenceville, GA 30043
USA
Phone: (877) 814-7004
Fax: (770) 234-3819
Web: www.vivebio.com
Sales: sales@vivebio.com
Support: support@vivebio.com

Authorized Representative:
AR-MED
A Division of Quintiles Ltd
Runnymede Malthouse Business Centre
off Hummer Road
Egham
Surrey TW20 9BD
United Kingdom

ANVISA: FRASCO E MATERIAL PARA COLETA ARMAZENAMENTO OU TRANSPORTE DE AMOSTRAS BIOLOGICAS
CAS 2535.247872/2010-95
CLASSE: I 805644740004
8436 CADASTRO DE PRODUTOS PARA DIAGNOSTICO DE USO IN-VITRO - IMPORTADO

ViveBio™, ViveST™ and A Hand in Life™ are registered trademarks of ViveBio, LLC. © 2011. All Rights Reserved.
No general patent(s) or license of any kind is granted other than this specific use from purchase of ViveST device.

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Test	AN concordantes	AN totales	Sensibilidad
TRUGENE	34.956,0	35.970,0	97%

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Referencias

1. RM Lloyd, Jr., DA Burns, JT Huong, RL Mathis, MA Winters, M Tanner, A De La Rosa, B Yen-Lieberman, W Armstrong, A Taege, DR McClernon, JL Wetshtein, BM Friedrich, MR Ferguson, W O'Brien, PM Feorino and M Holodniy. Dried-plasma transport using a novel matrix and collection system for human immunodeficiency virus and hepatitis C virus virologic testing. *J Clin Microbiol*, 47(5), May 2009, p. 1491-96.
2. M Zanoni, R Cortes, RS Diaz, MC Sucupira, D Ferreira, LA Inocencio, C Vilhena, C Loveday, RM Lloyd, Jr., M. Holodniy. Comparative effectiveness of dried plasma HIV-1 viral load testing in Brazil using ViveST for sample collection. *J Clin Virol*, 49(4), Dec 2010, p. 245-8.
3. National Committee for Clinical Laboratory Standards. Protection of laboratory workers from infectious disease transmitted by blood, body fluids, and tissue; approved guideline. NCCLS Document M29-A Villanova (PA); NCCLS; 1997 Dec. 90p.
4. Centers for Disease Control and Prevention: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other blood borne pathogens in healthcare settings. *MMWR*, 1988; 37: 377-82, 387-8.
5. United Nations: Recommendations on Transport of Dangerous Goods Model Regulations. UN document, 2009, section 2.6.3.
6. World Health Organization: Guidance on Regulations for the Transport of Infectious Substances 2009-2010. WHO document, 2008.

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Referências

1. RM Lloyd, Jr., DA Burns, JT Huang, RL Mathis, MA Winters, M Tanner, A De La Rosa, B Yen-Lieberman, W Armstrong, A Taege, DR McClermon, JL Wetshtein, BM Friedrich, MR Ferguson, W O'Brien, PM Feorino and M Holodniy. Dried-plasma transport using a novel matrix and collection system for human immunodeficiency virus and hepatitis C virus virologic testing. *J Clin Microbiol*, 47(5), May 2009, p. 1491-96.
2. M Zanoni, R Cortes, RS Diaz, MC Sucupira, D Ferreira, LA Inocencio, C Vilhena, C Loveday, RM Lloyd, Jr., M. Holodniy. Comparative effectiveness of dried plasma HIV-1 viral load testing in Brazil using ViveST for sample collection. *J Clin Virol*, 49(4), Dec 2010, p. 245-8.
3. National Committee for Clinical Laboratory Standards. Protection of laboratory workers from infectious disease transmitted by blood, body fluids, and tissue; approved guideline. NCCLS Document M29-A Villanova (PA); NCCLS; 1997 Dec. 90p.
4. Centers for Disease Control and Prevention: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other blood borne pathogens in healthcare settings. *MMWR*, 1988; 37: 377-82, 387-8.
5. United Nations: Recommendations on Transport of Dangerous Goods Model Regulations. UN document, 2009, section 2.6.3.
6. World Health Organization: Guidance on Regulations for the Transport of Infectious Substances 2009-2010. WHO document, 2008.

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ViveBio, LLC
1000 Hurricane Shoals Rd
Suite C330
Lawrenceville, GA 30043
EUA
Telefone: (877) 814-7004
Fax: (770) 234-3819
Página da Internet: www.vivebio.com

Representante autorizado:
AR-MED
A division of Quintiles Ltd
Runnymede Malthouse Business Centre
off Hummer Road
Egham
Surrey TW20 9BD
Reino Unido

Representante exclusivo no Brasil:
Centro de Genomas
Rua Afonso Celso, 469
CEP 04119-001
Sao Paulo/SP
Telefone: (11) 5079-9593
Página da Internet:
www.centrodegenomas.com.br

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ViveBio™, ViveST™ e A Hand in Life™ são marcas registadas da ViveBio, LLC. © 2011. Todos os direitos reservados.